


PM-LIFE® FIV IU

32

EMBRYO TRANSFER CATHETER



ECHOGENIC
Echogenic distal part for
better visualisation during
ultrasound



MANDREL
Used to pre-shape
the distal end

The PM-LIFE® FIV IU embryo transfer catheter is used for intrauterine insemination or to perform a micro-injection (transfer of embryos after in vitro fertilisation). Some product references come with an opaque tube for better ultrasound visualisation, or a mandrel to adapt the catheter to the patient's morphology.

ADVANTAGES

- Centimetric graduations and stopper for better control of catheter insertion.

AVAILABLE REFERENCE(S)

Reference	Catheter diameter (mm / Fr)	Inner diameter / volume (mm / ml)	Excess of the supple catheter (mm)	Useful length (mm)	Tube opaque	Mandrel
PS1280545	1.6 / 5	1.1 / 0.20	45	175	-	-
PS1280555	1.6 / 5	1.1 / 0.20	45	175	Yes	-
PS2380545	1.6 / 5	1.1 / 0.20	45	175	-	Yes

ABOUT THIS MEDICAL DEVICE

1 st release date	2002
Single use	Yes
Class	Is
Sterile	Yes
Sterilisation method	Sterilised via gamma radiation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 25 units in individual sterile wrapping
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
ZA La Sente Du Moulin
64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

PRINCE MEDICAL
DEVICES

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.